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**Metabolomics**

**Lipidomics**  
**Quantitation**

- Global untargeted Metabolomics
- Metabolite Quantitation

## Q Exactive Focus


**Quick and Easy**

- Detect
- Quantitate
- Confirm

# Q Exactive Focus: Great Reproducibility for Quan

D<sub>5</sub>-hippuric acid, external calibration, resolution = 72,000 at m/z 185

# Example 5: Rat Fasting Metabolism Study using FS-only with polarity switching

- Study designed to monitor the affect of fasting on metabolic profiles


Group	Male	Fasting Time <sup>1</sup>
1: 1101-1105	5	Dark Cycle Control (no Fast)
2. 2101-2105	5	2 hr Fast
3. 3101-3105	5	4 hr Fast
4. 4101-4105	5	8 hr Fast
5. 5101-5105	5	12 hr Fast
6. 6101-6105	5	16 hr Fast

<sup>1</sup> Rats fasted during a 6 p.m. to 6 a.m. dark cycle to capture peak feeding time

- **Samples:** 50uL Serum ppt with cold MeOH
- **MS:** 70K resolution, ESI+ and ESI-
- **Column:** Hypersil GOLD aQ 2.1x150mm, 1.9 $\mu$  @ 600 $\mu$ L/min, 50°C
- **Buffers:** A: 0.1% formic acid in H<sub>2</sub>O, B: 0.1% formic acid in 98:2 ACN:H<sub>2</sub>O


# Overall Method Robustness

## Uric Acid: Positive and Negative Data


# Translational Metabolomics Powered by TraceFinder 3.2 Software


- Productivity driven quantitation & screening
- SIEVE software plug-in with automatic import of compound DB
- Customized parameters for each metabolite
- Quantification of one compound based on calibration curve from another
- Efficient data review & powerful custom report generator

A selected compound viewed across all samples

All compounds viewed within a selected sample

Custom heat map by the new report builder plug-in

# Q Exactive Focus MS

## Challenge

Pharma

Biopharma

- Peptide Quan
- Biomarkers

## Q Exactive Focus

Quick and Easy

- Detect
- Quantitate
- Confirm


# Topic 6: HRAM Quantitation of Compounds in Plasma using targeted SIM analyses


- Readily Available Standards with labeled IS
  - Cerilliant, 1mg/mL in MeOH
- Samples spiked in rat plasma
  - Serial Dilution (1 pg/mL – 10,000pg/mL)
  - 4-5 Orders of magnitude linear dynamic range
  - Sample pooling (5 compounds in a single solution)
  - Replicate injections (n=6)
- Generic MS and LC conditions
  - MS - QE Focus – SIM @ 70K, generic source conditions
  - LC – U3000RS OAS, 8 min LC, flow 500uL/min, 5uL injection volume
  - Column – Hypersil Gold aQ, 2.1 x 50mm, 1.9um

# Verapamil Summary

Standard (pg/mL)	Mean Calc Amt	%CV	Count
Verapamil	--	--	0
5.00	--	--	0
10.0	10.4	11.1	4
50.0	44.2	5.1	6
100	88	1.3	6
500	478	2	6
1000	1011	1.7	6
5000	5199	0.9	6
10000	10820	0.9	6
50000	51944	2.3	6
100000	107655	3.1	6


RT: 4.15 - 4.66

NL:  
4.97E3  
m/z=  
455.28676-  
455.29404 MS  
ICIS  
Plasma\_Curve\_  
SIM\_MOD\_1\_3  
21

RT: 4.15 - 4.66

NL:  
4.97E3  
m/z=  
455.28676-  
455.29404 MS  
ICIS  
Plasma\_Curve\_  
SIM\_MOD\_1\_3  
21

# Summary Table

Compounds	Q Exactive Focus t-SIM (Res - 70,000)	
	LOD (ng/mL)	LOQ (ng/mL)
Buprenorphine	10	10
Paroxetine	5	10
Verapamil	10	10
Alprazolam-D5	10	10
Reserpine	5	50

# Conclusion

- Main requirements for quantitation
  - Sensitivity ✓
  - Linear Response ✓
  - Reproducibility ✓
- Ruggedness/Robustness
  - Preventative Maintenance ✓
  - Method Transfer ✓
  - Flexibility ✓
- Method Development/Troubleshooting
  - Orthogonal Technique ✓
  - Little to no optimization needed ✓
  - Interference and background characterization ✓

# Q Exactive Focus Applications Universe

# Q Exactive Focus Software Solutions

## TraceFinder™ 3.2

Optimized for Forensic Toxicology

© Copyright 2013 Thermo Fisher Scientific Inc.  
All rights reserved. This program is protected by copyright  
law and international treaties as described in Help About.

## TraceFinder™ 3.2

Optimized for Environmental and Food Safety

© Copyright 2013 Thermo Fisher Scientific Inc.  
All rights reserved. This program is protected by copyright  
law and international treaties as described in Help About.

Thermo ToxFinder 1.0

## ToxFinder 1.0

Full Scan -- AIF

Full Scan -- Data Dependent

SRM

© Copyright 2014 Thermo Fisher Scientific Inc. All rights reserved. This program is protected by copyright law and international treaties as described in Help About.

**Thermo**  
SCIENTIFIC

# PORTFOLIO

## Exactive Plus (EMR)

- Orbitrap analyzer
- Mass Range m/z 50 - 6000
- Mass Range EMR: 300 - 20,000
- Mass Accuracy: <1ppm
- Mass Resolution >140,000
- Scan Speed up to 12 Hz

## Q Exactive Focus

- Orbitrap analyzer
- Mass Range m/z 50 - 2000
- Mass Accuracy <1ppm
- Max. Mass Resolution >70,000
- Scan speed up to 12Hz MS and 12 Hz MS/MS
- Polarity switching
- Top 2 ddMS2
- No Multiplexing
- Refined workflows that are optimal for routine laboratories

## Q Exactive & Q Exactive Plus

- Orbitrap analyzer
- Mass Range m/z 50 - 6000
- Mass Accuracy <1ppm
- Max. Mass Resolution >140,000
- Scan speed up to 12Hz MS and MS/MS
- Spectral Multiplexing
- AQT & AABG (QE Plus only)
- Optional Intact Protein Mode and Enhanced Resolution (280k) for QE Plus only

## Q Exactive HF

- Ultra High Field Orbitrap analyzer
- Mass Range m/z 50 - 6000
- Mass Accuracy <1ppm
- Max. Mass Resolution >240,000
- Scan speed up to 18Hz MS and MS/MS
- Spectral Multiplexing
- AQT & AABG
- Optional Protein Mode

# Overview of Current Orbitrap MS Portfolio

							
Attributes	Exactive Plus (EMR)	Q Exactive Focus	Q Exactive (Plus)	Q Exactive HF	Orbitrap LTQ XL	Orbitrap Elite	Orbitrap Fusion
Max Resolution (FWHM) @m/z 200	140k	70k	140k (280k Enh R)	240k	140k	340k	500k
Mass accuracy, ppm (internal cal)	<1	<1	<1	<1	<1	<1	<1
Scan rate (Hz)	12 Hz	12Hz	12 Hz	18Hz	1 Hz	4 Hz	18 Hz
Polarity switch (s)	<1	<1	<1	<1	Not during run	Not during run	1.1
Parallel Reaction Monitoring	NA	Yes	Yes	Yes	Yes	Yes	Yes
Multiplex (precursor/scan)	NA	No	10	10	NA	NA	10
DIA-MS / SWATH	NA	vDIA	SWATH™-like, vDIA, msxDIA	SWATH™-like, vDIA, msxDIA	NA	Same as OT Fusion	WiSIM-DIA, SWATH™-like, vDIA, msxDIA
MS <sup>n</sup> (n=10)	NA	NA	NA	NA	Yes	Yes	Yes
ETD option	NA	NA	NA	NA	Yes	Yes	Yes
Decision-tree (CID/HCD/ETD)	NA	NA	NA	NA	NA	Yes	Yes
Synchronous MS <sup>3</sup>	NA	NA	NA	NA	NA	NA	Yes

**CONCLUSION**

# Top Five Reasons to Buy

1

**Higher performance analysis in the most affordable Orbitrap ever, refined for routine laboratory experiments**

2

**Use for routine Qual and Quan with optional retrospective analysis, two techniques on one platform**

3

**Easy to use; quick method set up, no optimization**

4

**Versatile Processing Software for both Qual & Quan workflows**

5

**User Friendly in Daily Operation**

# Marketing Collateral: Q Exactive Focus MS

- Brochure
- Specification Sheet
- Application and Technical Notes
- Animation
- Testimonials
- Scientific Presentations
- Planet Orbitrap and RevBase promotions

**Detection of Stanozolol Glucuronides in Human Sports Drug Testing by Means of High-Resolution, Accurate-Mass Mass Spectrometry**

**Quantitation of Opiates to Low ng/mL Levels in Urine for Forensic Use Using an Affordable, High-Resolution, Accurate-Mass Mass Spectrometer**

**Metabolomic Profiling in Drug Discovery: Understanding the Factors that Influence a Metabolomics Study and Strategies to Reduce Biochemical and Chemical Noise**

**Quick and Sensitive Analysis of Multiclass Veterinary Drug Residues in Meat, Plasma, and Milk on a Q Exactive Focus LC-MS System**

**Thermo Scientific Q Exactive Focus Orbitrap LC-MS/MS System**  
Affordable. Durable. Proven.

- Scan speed up to 12 Hz with Orbitrap Analyser Technology for best sensitivity and quantitative results
- Resolving power of up to 70,000 (FWHM) at m/z 200
- Routine sub ppm mass accuracy
- Linear Dynamic Range up to 6 orders of magnitude
- Multiple approaches to quantitation including Selected Ion Monitoring (SIM), Parallel Reaction Monitoring (PRM), and Data Independent Acquisition (DIA)
- Priority scheduling for maximum instrument coverage
- High-Energy Collisional Dissociation (HCD)

**The Thermo Scientific™ Q Exactive™ Focus Orbitrap LC-MS/MS combines quadrupole precursor ion selection and a high-resolution accurate mass Orbitrap Orbitrap mass analyser to deliver mass accuracy, sensitivity, and fast polarity switching while maintaining usability comparable to a triple quadrupole mass spectrometer.**

**The Q Exactive Focus mass spectrometer delivers qualitative and quantitative results every day. Perform routine identification, quantitation and confirmation in a single analysis with easy-to-use Thermo Scientific™ TracePro™ acquisition and processing software.**

**Based on the proven Q Exactive platform, the Q Exactive Focus provides durable and reliable performance for routine laboratory drug and metabolite screening, profiling and quantitative analysis of small molecule compounds. Full scan confirmation mode and parallel reaction monitoring (PRM) provide reproducible quantitation results and targeted sensitivity capabilities. High-resolution independent analysis (DIA) provides complete qualitative coverage for unbiased screening without compromising precise quantitative attributes.**

**The Q Exactive Focus system provides best-in-class specificity at value price. It's ideal for routine labs performing food safety residue analysis, environmental analysis, forensic toxicology, sports doping, clinical research, metabolomics, and pharmaceutical research.**

**Thermo Scientific**

The image features a dark blue background with vertical light streaks and a horizontal light flare. A white sine wave with small starburst highlights is positioned below the company name.

**Thermo**  
S C I E N T I F I C

**Transform Your Science**